


1653 East Main Street  
Rochester, NY 14609 USA  
Voice: 585.482.0300  
FAX: 585.288.5989  
imaging@appliedimage.com

T-22  
Improved Layout  
USAF 1951 Chart  
Product Spec Sheet




Catalog Part No: (grouped by substrate material)

**T-22-P-CG / T-22-2-N-CG**  
**T-22-P-OP / T-22-2-N-OP**  
**T-22-N-TM / T-22-P-TM**  
**T-22-P-RM**

Product Name: **USAF 1951 Chart**, (improved layout and features)


Drawing / Photo of Part:



**This chart appears to be different than the USAF-1951 Standard Chart but is easier to use and meets all requirements of MIL-STD-150A!**

The above image is an approximate representation of the actual product.  
Specifications are subject to change without notice.

Description: This test target meets requirements as specified in MIL-STD-150A for resolving power tests since contrast and feature sizes are its only requirements. All bars and space features are the same as provided on the traditional USAF-1951 Chart (T-20 USAF-1951 Chart). Target *Elements* are arranged in order of increasing frequency from top to bottom in two columns, with frequencies of 0.25 cycles per mm (c/mm) through 0.45 c/mm in the left column and frequencies of 0.50 c/mm and higher in the right column. In order to avoid the necessity of a lookup chart and for maximum recognizability, a label using an OCR font to indicate frequency in c/mm, is adjacent to each *Element*. An 8x8mm square is provided to aid in density measurements. This version of the USAF standard chart

1653 East Main Street Rochester, NY 14609 USA Voice: 585.482.0300 FAX: 585.288.5989 imaging@appliedimage.com	T-22 Improved Layout USAF 1951 Chart Product Spec Sheet	APPLIED® IMAGE Inc 
--	--	--

has *Finder Square* features to better identify elements when the chart is used near the limit of resolution. *Finder Squares* are provided in the following table. (Note: Applied Image offers the standard USAF-1951 chart as part number T-20.)

<i>Element Frequency</i>	<i>Number of Squares</i>	<i>Element Frequency (c/mm)</i>	<i>Number of Squares</i>
1.0 (c/mm)	1	20.1 (c/mm)	2
2.0 (c/mm)	2	50.8 (c/mm)	5
5.0 (c/mm)	5	101 (c/mm)	1
10.1 (c/mm)	1	203 (c/mm)	2

**Substrate Size, Type, Image Forming Material, Polarity and Frequency Range:**

<i>Part Number</i>	<i>Substrate Size (mm)</i>	<i>Substrate Type</i>	<i>Image Forming Material</i>	<i>Polarity</i>	<i>Frequencies (c/mm)</i>
T-22-P-CG	102 x 102	Soda-lime glass	Chrome	Positive	0.25 → 228
T-22-P-OP	102 x 102	White opal glass	Chrome	Positive	0.25 → 228
T-22-2-N-CG	102 x 102	Soda-lime glass	Chrome	Negative	0.25 → 512
T-22-2-N-OP	102 x 102	White opal glass	Chrome	Negative	0.25 → 512
T-22-N-TM	102 x 102	Clear polyester	Photo-emulsion	Negative	0.25 → 128
T-22-P-TM	102 x 102	Clear polyester	Photo-emulsion	Positive	0.25 → 128
T-22-P-RM	102 x 102	Photo-paper	Photo-emulsion	Positive	0.25 → 22.6

Reading Direction: Right Read Chrome / Emulsion Up (RRCU / RREU).

Image Placement Accuracy: 0.001mm/100mm on glass. Polyester or paper expansion is the major source of change in photo-emulsion versions.

Feature Size Accuracy: Typical line width is 2 to 5 percent of aim.

Image Contrast / Density: Chrome transmission density is 2.0 or higher; photo-emulsion transmission density 2.0 or higher (photo-emulsion on polyester); emulsion reflection density 1.3 or higher (photo-emulsion on paper).

History / typical use: Checking of optical resolving power for optical systems.

How To Use Resolution Test Charts: Multi-bar test charts are most often used to determine the overall resolution limit of an entire optical system. Each component of a system adds some degradation to image quality. Therefore, the comparison of a test chart in the object plane to the resulting image (image plane vs. film or image vs. electronic file) will produce the complete account of all components such as optics, image processing, vibration, color fringing, etc. Test charts may be used in optical systems that reduce, maintain (the same image size 1:1) or magnify

1653 East Main Street  
Rochester, NY 14609 USA  
Voice: 585.482.0300  
FAX: 585.288.5989  
imaging@appliedimage.com

T-22  
Improved Layout  
USAF 1951 Chart  
Product Spec Sheet




the image size however; system magnification must be applied as a multiplier to determine the resolving power limit. **Example:** To evaluate a 35mm film-camera system, load the camera with the film of interest. Position the camera at a measured distance from the test chart and achieve good focus. (Commonly 25 focal lengths is the distance used from camera to the test chart.) Develop the image and measure the object and image sizes. (i.e. If the test chart card is 100mm wide in the object plane and 5.5mm on the film, the system reduction factor is 1:18.18 (0.055). Next, examine the film image using a microscope and determine the resolution limit. If the last distinct *Element* has a frequency of 5.04 c/mm, then the resolving power or the resolution of the film is calculated from that image frequency divided by the reduction (5.04/0.055 or 91.6 c/mm).

How to read a resolution bar test pattern: Resolution is defined as a threshold measurement. The resolution limit is the last group of bars that can be seen and correctly counted in both horizontal and vertical directions. Breaks within the bar pattern are acceptable, and do not cause failure of resolution as long as the bars are distinctly identified and accurately counted. A microscope with the necessary magnification is commonly used to view reduced images. The correct number of bars must be counted when determining the resolution limit in order to avoid an inflated measure. This problem can be seen while observing increasingly finer groups. Features will be beyond resolution and then appear again. One cause of this spurious resolution is poor focus of a lens, which causes the bars and spaces in the image to apparently flip polarity. Other causes are artifacts from digital scanning, or image processing algorithms.

Target Viewing Angle: Test charts may be placed at the center and outer edges of the field of view so that the comparisons can be made of resolving power between locations across the field of view. Commonly, a comparison is made between center and edge at 90-degree increments at the outer edge of the field of view. This is useful for detection of directional vibration influences and the common problem of uniformity of quality throughout the lens field of view. A more complete assessment would use charts at additional locations throughout the field for example, comparing resolving power at the center with the entire outside edge at 45-degree increments.

1653 East Main Street  
Rochester, NY 14609 USA  
Voice: 585.482.0300  
FAX: 585.288.5989  
imaging@appliedimage.com

T-22  
Improved Layout  
USAF 1951 Chart  
Product Spec Sheet



Terms: Line Pairs per mm (LP/mm) or lines per mm (lpmm) equal cycles per mm.


Related documents:

c-mm LPI DPI info.pdf  
T-20\_spec.pdf  
T-21\_spec.pdf

Related Parts: T-21, T-20, T-10, T-90

1653 East Main Street  
 Rochester, NY 14609 USA  
 Voice: 585.482.0300  
 FAX: 585.288.5989  
 imaging@appliedimage.com


T-22  
 Improved Layout  
 USAF 1951 Chart  
 Product Spec Sheet



Group	Element	Frequency (cycles/mm)	Line or space width (micron)	Line or space width (inches)	DPI (Dots per inch)	PPI (Pixels per Inch)	LPPI (line pairs per inch)	Dots per Cm	Lines pairs per Cm
-6	1	0.016	32000	1.259843		0.79	0.40	0.31	0.16
	2	0.018	28509	1.122392		0.89	0.45	0.35	0.18
	3	0.020	25398	0.999938		1.00	0.50	0.39	0.20
	4	0.022	22627	0.890843		1.12	0.56	0.44	0.22
	5	0.025	20159	0.793651		1.26	0.63	0.50	0.25
	6	0.028	17959	0.707063		1.41	0.71	0.56	0.28
-5	1	0.031	16000	0.629921		1.59	0.79	0.63	0.31
	2	0.035	14254	0.561196		1.78	0.89	0.70	0.35
	3	0.039	12699	0.499969		2.00	1.00	0.79	0.39
	4	0.044	11314	0.445422		2.25	1.12	0.88	0.44
	5	0.050	10079	0.396826		2.52	1.26	0.99	0.50
	6	0.056	8980	0.353531		2.83	1.41	1.11	0.56
-4	1	0.063	8000	0.314961		3.18	1.59	1.25	0.63
	2	0.070	7127	0.280598		3.56	1.78	1.40	0.70
	3	0.079	6350	0.249984		4.00	2.00	1.57	0.79
	4	0.088	5657	0.222711		4.49	2.25	1.77	0.88
	5	0.099	5040	0.198413		5.04	2.52	1.98	0.99
	6	0.111	4490	0.176766		5.66	2.83	2.23	1.11
-3	1	0.125	4000	0.157480		6.35	3.18	2.50	1.25
	2	0.140	3564	0.140299		7.13	3.56	2.81	1.40
	3	0.157	3175	0.124992		8.00	4.00	3.15	1.57
	4	0.177	2828	0.111355		8.98	4.49	3.54	1.77

1653 East Main Street  
 Rochester, NY 14609 USA  
 Voice: 585.482.0300  
 FAX: 585.288.5989  
 imaging@appliedimage.com

T-22  
 Improved Layout  
 USAF 1951 Chart  
 Product Spec Sheet




Group	Element	Frequency (cycles/mm)	Line or space width (micron)	Line or space width (inches)	DPI (Dots per inch)	PPI (Pixels per Inch)	LPPI (line pairs per inch)	Dots per Cm	Lines pairs per Cm
	5	0.198	2520	0.099206		10.08	5.04	3.97	1.98
	6	0.223	2245	0.088383		11.31	5.66	4.45	2.23
-2	1	0.250	2000	0.078740		12.70	6.35	5.00	2.50
	2	0.281	1782	0.070150		14.26	7.13	5.61	2.81
	3	0.315	1587	0.062496		16.00	8.00	6.30	3.15
	4	0.354	1414	0.055678		17.96	8.98	7.07	3.54
	5	0.397	1260	0.049603		20.16	10.08	7.94	3.97
	6	0.445	1122	0.044191		22.63	11.31	8.91	4.45
-1	1	0.500	1000.00	0.039370		25.40	12.70	10.00	5.00
	2	0.561	890.90	0.035075		28.51	14.26	11.22	5.61
	3	0.630	793.70	0.031248		32.00	16.00	12.60	6.30
	4	0.707	707.11	0.027839		35.92	17.96	14.14	7.07
	5	0.794	629.96	0.024802		40.32	20.16	15.87	7.94
	6	0.891	561.23	0.022096		45.26	22.63	17.82	8.91
0	1	1.000	500.00	0.019685		50.80	25.40	20.00	10.00
	2	1.122	445.45	0.017537		57.02	28.51	22.45	11.22
	3	1.260	396.85	0.015624		64.00	32.00	25.20	12.60
	4	1.414	353.55	0.013919		71.84	35.92	28.28	14.14
	5	1.587	314.98	0.012401		80.64	40.32	31.75	15.87
	6	1.782	280.62	0.011048		90.52	45.26	35.64	17.82
1	1	2.000	250.00	0.009843		101.60	50.80	40.00	20.00
	2	2.245	222.72	0.008769		114.04	57.02	44.90	22.45

Please contact Applied Image customer service at the address noted above, for custom images, shapes and materials.

1653 East Main Street  
 Rochester, NY 14609 USA  
 Voice: 585.482.0300  
 FAX: 585.288.5989  
 imaging@appliedimage.com

T-22  
 Improved Layout  
 USAF 1951 Chart  
 Product Spec Sheet




Group	Element	Frequency (cycles/mm)	Line or space width (micron)	Line or space width (inches)	DPI (Dots per inch)	PPI (Pixels per Inch)	LPPI (line pairs per inch)	Dots per Cm	Lines pairs per Cm
	3	2.520	198.43	0.007812	128.01	64.00	50.40	25.20	
	4	2.828	176.78	0.006960	143.68	71.84	56.57	28.28	
	5	3.175	157.49	0.006200	161.28	80.64	63.50	31.75	
	6	3.564	140.31	0.005524	181.03	90.52	71.27	35.64	
2	1	4.00	125.00	0.004921	203.20	101.60	80.00	40.00	
	2	4.49	111.36	0.004384	228.08	114.04	89.80	44.90	
	3	5.04	99.21	0.003906	256.02	128.01	100.79	50.40	
	4	5.66	88.39	0.003480	287.37	143.68	113.14	56.57	
	5	6.35	78.75	0.003100	322.56	161.28	126.99	63.50	
	6	7.13	70.15	0.002762	362.06	181.03	142.54	71.27	
3	1	8.00	62.50	0.002461	406.40	203.20	160.00	80.00	
	2	8.98	55.68	0.002192	456.17	228.08	179.59	89.80	
	3	10.08	49.61	0.001953	512.03	256.02	201.59	100.79	
	4	11.31	44.19	0.001740	574.74	287.37	226.27	113.14	
	5	12.70	39.37	0.001550	645.12	322.56	253.98	126.99	
	6	14.25	35.08	0.001381	724.12	362.06	285.09	142.54	
4	1	16.00	31.25	0.001230	812.80	406.40	320.00	160.00	
	2	17.96	27.84	0.001096	912.34	456.17	359.19	179.59	
	3	20.16	24.80	0.000977	1024.06	512.03	403.17	201.59	
	4	22.63	22.10	0.000870	1149.47	574.74	452.55	226.27	
	5	25.40	19.69	0.000775	1290.24	645.12	507.97	253.98	
	6	28.51	17.54	0.000690	1448.24	724.12	570.18	285.09	

Please contact Applied Image customer service at the address noted above, for custom images, shapes and materials.

1653 East Main Street  
 Rochester, NY 14609 USA  
 Voice: 585.482.0300  
 FAX: 585.288.5989  
 imaging@appliedimage.com

T-22  
 Improved Layout  
 USAF 1951 Chart  
 Product Spec Sheet




Group	Element	Frequency (cycles/mm)	Line or space width (micron)	Line or space width (inches)	DPI (Dots per inch)	PPI (Pixels per Inch)	LPPI (line pairs per inch)	Dots per Cm	Lines pairs per Cm
5	1	32.00	15.63	0.000615	1625.60	812.80		640.00	320.00
	2	35.92	13.92	0.000548	1824.67	912.34		718.38	359.19
	3	40.32	12.40	0.000488	2048.13	1024.06		806.35	403.17
	4	45.25	11.05	0.000435	2298.95	1149.47		905.10	452.55
	5	50.80	9.84	0.000388	2580.48	1290.24		1015.94	507.97
	6	57.02	8.77	0.000345	2896.49	1448.24		1140.35	570.18
6	1	64.00	7.81	0.000308	3251.20	1625.60		1280.00	640.00
	2	71.84	6.96	0.000274	3649.35	1824.67		1436.75	718.38
	3	80.63	6.20	0.000244	4096.26	2048.13		1612.70	806.35
	4	90.51	5.52	0.000217	4597.89	2298.95		1810.19	905.10
	5	101.59	4.92	0.000194	5160.96	2580.48		2031.87	1015.94
	6	114.04	4.38	0.000173	5792.98	2896.49		2280.70	1140.35
7	1	128.00	3.91	0.000154	6502.40	3251.20		2560.00	1280.00
	2	143.68	3.48	0.000137	7298.70	3649.35		2873.50	1436.75
	3	161.27	3.10	0.000122	8192.51	4096.26		3225.40	1612.70
	4	181.02	2.76	0.000109	9195.78	4597.89		3620.39	1810.19
	5	203.19	2.46	0.000097	10321.92	5160.96		4063.75	2031.87
	6	228.07	2.19	0.000086	11585.96	5792.98		4561.40	2280.70
8	1	256.00	1.95	0.000077	13004.80	6502.40		5120.00	2560.00
	2	287.35	1.74	0.000069	14597.39	7298.70		5747.01	2873.50
	3	322.54	1.55	0.000061	16385.02	8192.51		6450.80	3225.40
	4	362.04	1.38	0.000054	18391.56	9195.78		7240.77	3620.39



1653 East Main Street  
 Rochester, NY 14609 USA  
 Voice: 585.482.0300  
 FAX: 585.288.5989  
 imaging@appliedimage.com

T-22  
 Improved Layout  
 USAF 1951 Chart  
 Product Spec Sheet



Group	Element	Frequency (cycles/mm)	Line or space width (micron)	Line or space width (inches)	DPI (Dots per inch)	PPI (Pixels per Inch)	LPPI (line pairs per inch)	Dots per Cm	Lines pairs per Cm
	5	406.37	1.23	0.000048	20643.83	10321.92		8127.49	4063.75
	6	456.14	1.10	0.000043	23171.92	11585.96		9122.80	4561.40
9	1	512.00	0.98	0.000038	26009.60	13004.80		10240.00	5120.00
	2	574.70	0.87	0.000034	29194.79	14597.39		11494.01	5747.01
	3	645.08	0.78	0.000031	32770.04	16385.02		12901.59	6450.80
	4	724.08	0.69	0.000027	36783.13	18391.56		14481.55	7240.77
	5	812.75	0.62	0.000024	41287.67	20643.83		16254.99	8127.49
	6	912.28	0.55	0.000022	46343.84	23171.92		18245.61	9122.80